

TransIT®-Jurkat Transfection Reagent

Quick Reference Protocol

Instructions for MIR 2120, 2122, 2124, 2125, 2126

Full protocol, SDS and Certificate of Analysis available at mirusbio.com/2120



SPECIFICATIONS

Storage	Store TransIT®-Jurkat Reagent tightly capped at 4°C. <i>Before each use</i> , warm to room temperature and vortex gently.
Product Guarantee	1 year from the date of purchase, when properly stored and handled.

► PLASMID DNA TRANSFECTION PROTOCOL



Full protocol and additional documentation available at mirusbio.com/2120

Fill in volumes below based on culture vessel used for transfection (Table 1).

- A. Plate cells**
- 1. Plate cells in ____ ml complete growth medium (per well) at a density of $2-4 \times 10^5$ cells/ml.
 - 2. Culture overnight.
Optional: Alternatively, plate cells at a density of $4-8 \times 10^5$ cells/ml complete growth medium, just prior to transfection.
- B. Prepare TransIT®-Jurkat Reagent:DNA complexes**
- 1. Warm TransIT®-Jurkat to room temperature and vortex gently.
 - 2. Place ____ µl of Opti-MEM® I Reduced-Serum Medium in a sterile tube.
 - 3. Add ____ µl plasmid DNA. Mix gently by pipetting.
 - 4. Add ____ µl of TransIT®-Jurkat Reagent. Mix gently by pipetting.
 - 5. Incubate at room temperature for 15-30 minutes.
- C. Distribute complexes to cells**
- 1. Add TransIT®-Jurkat:DNA complex mixture drop-wise to different areas of the well.
 - 2. Gently rock plate for even distribution of complexes.
 - 3. Incubate 24-72 hours.
 - 4. Harvest cells and assay as required.

Table 1. Recommended starting conditions

Culture vessel	24-well plate	12-well plate	6-well plate
Surface area	1.9 cm ²	3.8 cm ²	9.6 cm ²
Complete growth medium	0.5 ml	1 ml	2.5 ml
Serum-free medium	50 µl	100 µl	250 µl
DNA (1 µg/µl stock)	0.5 µl	1 µl	2.5 µl
TransIT®-Jurkat Reagent	1.5 µl	3 µl	7.5 µl

► Transfection Optimization

Determine the best TransIT®-Jurkat Reagent:DNA ratio for each cell type. Start with 3 µl of TransIT®-Jurkat Reagent per 1 µg of DNA. Vary the concentration of TransIT®-Jurkat Reagent from 1-5 µl per 1 µg DNA to find the optimal ratio.

TransIT®-Jurkat reagent also works well for additional cell lines of hematopoietic origin such as K562, RAW264.7 and THP-1.

For additional optimization tips, see [full protocol](#).

► NOTES



Reagent Agent®

Reagent Agent® is an online tool designed to help determine the best solution for nucleic acid delivery based on in-house data, customer feedback and citations.

Learn more at: mirusbio.com/ra

©1996-2025 All rights reserved. Mirus Bio LLC. All trademarks are the property of their respective owners.

For terms and conditions, visit www.mirusbio.com

Rev.C 042721

Mirus Bio LLC

www.mirusbio.com | techsupport@mirusbio.com | U.S. Toll Free: 844.MIRUSBIO (844.647.8724) | Direct: +1.608.441.2852