

TransIT-PRO® Transfection Reagent for Protein
Production in Expi293F™ Cells



Quick Reference Protocol
Instructions for MIR 5720, 5730, 5740 and 5750

SPECIFICATIONS

Storage	Store TransIT-PRO® Reagent tightly capped at -20°C. <i>Before each use</i> , warm to room temperature and vortex gently.
Product Guarantee	1 year from the date of purchase, when properly stored and handled.

► PLASMID DNA TRANSFECTION PROTOCOL



Full protocol and additional documentation available at mirusbio.com/5740

Fill in volumes below based on total culture volume (Table 1).

A. Maintenance of cells

- 1. Passage Expi293F™ cells to a density of 2×10^6 cells/ml in Expi293™ Expression Medium 18–24 hours prior to transfection to ensure that cells are actively dividing at the time of transfection.
- 2. Culture cells overnight on an orbital shaker (125 rpm when using a shaker with a 2 cm orbital throw) at appropriate temperature and CO₂ levels (e.g. 37°C, 5-8% CO₂).

B. Prepare TransIT-PRO® Reagent:DNA complexes

NOTE: This procedure recommends a 1:1 TransIT-PRO® Reagent:DNA ratio, which is optimal for most proteins. However, doubling both DNA and reagent amounts (i.e. 60 µl TransIT-PRO® Reagent + 60 µg total DNA per 30 ml culture volume) may result in improved expression with some constructs.

- 1. Determine cell density and viability. DO NOT proceed with transfection if cells are <95% viable.
- 2. Dilute Expi293F™ cells to a final density of 2.5×10^6 cells/ml in Expi293™ Medium immediately before transfection. Return cells to the incubator while preparing transfection complexes.
- 3. Warm TransIT-PRO® Reagent to room temperature and vortex gently.
- 4. Place ____ ml of serum-free medium (e.g. Opti-MEM® or Opti-PRO™ SFM) in a sterile tube.
- 5. Add ____ µg plasmid DNA. Mix gently by pipetting.
- 6. Add ____ µl of TransIT-PRO® Reagent. Mix gently by pipetting.
- 7. Incubate complexes at room temperature for 15-20 minutes.

C. Distribute complexes to cells

- 1. Add TransIT-PRO® Reagent:DNA complexes to cultured Expi293F™ cells.
- 2. Culture transfected cells at appropriate conditions (e.g. 37°C, 5-8% CO₂, shaking) until harvest. The optimal incubation time for protein expression and harvest depends on the expression construct, nature of the protein, and detection method (i.e. 2-14 days post-transfection).
- 3. Harvest cells and/or supernatant and assay as required.

Table 1. Volume scaling worksheet for DNA transfections with TransIT-PRO® Transfection Reagent.

Starting conditions per milliliter of complete growth medium				
	Per 1 ml		Total culture volume	Reagent quantities
Serum-free Complex Medium	0.1 ml	×	_____ ml	= _____ ml
Plasmid DNA (1 µg/µl stock)	1 µl	×	_____ ml	= _____ µl
TransIT-PRO® Transfection Reagent	1 µl	×	_____ ml	= _____ µl

► Additional Information

Expi293™ and Expi293F™ are trademarks of Life Technologies Corporation.
For additional optimization tips, see TransIT-PRO® Transfection Reagent [full protocol](#).

Mirus Bio LLC

www.mirusbio.com | techsupport@mirusbio.com | U.S. Toll Free: 844.MIRUSBIO (844.647.8724) | Direct: +1.608.441.2852



Reagent Agent[®]

Reagent Agent[®] is an online tool designed to help determine the best solution for nucleic acid delivery based on in-house data, customer feedback and citations.

Learn more at: mirusbio.com/ra

©1996-2024 All rights reserved. Mirus Bio LLC. All trademarks are the property of their respective owners.
For terms and conditions, visit www.mirusbio.com

Rev.A 042920