

TransIT-TKO® Transfection Reagent

Quick Reference Protocol

Instructions for MIR 2150, 2152, 2154, 2155, 2156
Full protocol, SDS and Certificate of Analysis available at mirusbio.com/2150



SPECIFICATIONS

Storage	Store TransIT-TKO® Reagent tightly capped at 4°C. <i>Before each use</i> , warm to room temperature and vortex gently.
Product Guarantee	1 year from the date of purchase, when properly stored and handled.

► siRNA TRANSFECTION PROTOCOL



Full protocol and additional documentation available at mirusbio.com/2150

Fill in volumes below based on culture vessel used for transfection (Table 1).

- A. Plate cells
- 1. Plate cells in ___ml complete growth medium (per well).
For adherent cells: Plate cells at a density of $0.8-3.0 \times 10^5$ cells/ml.
For suspension cells: Plate cells at a density of $2.5-5.0 \times 10^5$ cells/ml.
 - 2. Culture overnight. Most cell types should be $\geq 80\%$ confluent on day of transfection.
- B. Prepare TransIT-TKO® Reagent:siRNA complexes
- 1. Warm TransIT-TKO® to room temperature and vortex gently.
 - 2. Place ___µl of Opti-MEM® I Reduced-Serum Medium in a sterile tube.
 - 3. Add ___µl TransIT-TKO® Reagent. Mix gently by pipetting.
 - 4. Add ___µl of a 10 µM siRNA stock solution (25 nM final concentration). Mix gently by pipetting.
 - 5. Incubate at room temperature for 15-30 minutes.
- C. Distribute complexes to cells
- 1. Add TransIT-TKO® Reagent:siRNA complex mixture drop-wise to different areas of the well.
 - 2. Gently rock plate for even distribution of complexes.
 - 3. Incubate 24-72 hours.
 - 4. Harvest cells and assay for knockdown of gene expression.

Table 1. Recommended starting conditions

Culture vessel	24-well plate	12-well plate	6-well plate
Surface area	1.9 cm ²	3.8 cm ²	9.6 cm ²
Complete growth medium	0.5 ml	1 ml	2.5 ml
Serum-free medium	50 µl	100 µl	250 µl
TransIT-TKO® Reagent	2.5 µl	5 µl	10 µl
siRNA (10 µM stock, 25 nM final)	1.4 µl	2.8 µl	6.8 µl

► Transfection Optimization

Determine the best volume of TransIT-TKO® for each cell type. Start with 10 µl of TransIT-TKO® per well of a 6-well plate. For further optimization, vary the amount from 8-12 µl per well to find the optimal volume.

For more tips and instructions for co-transfection, see full protocol. Cell-type-specific recommendations available at: Reagent Agent: mirusbio.com/ra



Reagent Agent®

Reagent Agent® is an online tool designed to help determine the best solution for nucleic acid delivery based on in-house data, customer feedback and citations.

Learn more at: mirusbio.com/ra

©1996-2025 All rights reserved. Mirus Bio LLC. All trademarks are the property of their respective owners.

For terms and conditions, visit www.mirusbio.com

Rev.C 042721