

# TransIT-X2® Dynamic Delivery System

## Quick Reference Protocol

[DNA Delivery Instructions for MIR 6000, 6003, 6004, 6005, 6006, 6010](#)

Full protocol, SDS and Certificate of Analysis available at [mirusbio.com/6000](http://mirusbio.com/6000)



## SPECIFICATIONS

Storage	Store TransIT-X2® Dynamic Delivery System tightly capped at -20°C. <b>Before each use</b> , warm to room temperature and vortex gently.
Product Guarantee	1 year from the date of purchase, when properly stored and handled.

### ▶ TRANSFECTION PROTOCOL | PLASMID DNA



Full protocol and additional documentation available at [mirusbio.com/6000](http://mirusbio.com/6000)

### Fill in volumes below based on culture vessel used for transfection (Table 1).

#### A. Plate cells

1. Plate cells in \_\_\_ ml complete growth medium (per well or flask) 18-24 hours prior to transfection to ensure cells are actively dividing at the time of transfection.

**For adherent cells:** Plate cells at a density of 0.8—3.0 x 10<sup>5</sup> cells/ml.

**For suspension cells:** Plate cells at a density of 2.5—5.0 x 10<sup>5</sup> cells/ml.

2. Culture overnight. Most cell types should be ≥80% confluent at time of transfection.

#### B. Prepare TransIT-X2® Reagent:DNA complexes

1. Warm TransIT-X2® Reagent to room temperature and vortex gently.
2. Place \_\_\_ μl of OptiMEM® I Reduced-Serum Medium in a sterile tube.
3. Add \_\_\_ μl plasmid DNA. Mix gently by pipetting.
4. Add \_\_\_ μl of TransIT-X2®. Mix gently by pipetting.
5. Incubate at room temperature for 15-30 minutes.

#### C. Distribute complexes to cells

1. Add TransIT-X2®:DNA complex mixture drop-wise to different areas of the well.
2. Gently rock plate for even distribution of complexes.
3. Incubate 24-72 hours.
4. Harvest cells and assay as required.

Table 1. Recommended starting conditions

Culture vessel	24-well plate	12-well plate	6-well plate
Surface area	1.9 cm <sup>2</sup>	3.8 cm <sup>2</sup>	9.6 cm <sup>2</sup>
Complete growth medium	0.5 ml	1 ml	2.5 ml
Serum-free medium	50 μl	100 μl	250 μl
DNA (1 μg/μl stock)	0.5 μl	1 μl	2.5 μl
TransIT-X2® Reagent	1.5 μl	3 μl	7.5 μl

### ▶ Transfection Optimization

Determine the best TransIT-X2®:DNA ratio for each cell type. Start with 3 μl of TransIT-X2® per 1 μg of DNA. Vary the concentration of TransIT-X2® from 2–6 μl per 1 μg DNA to find the optimal ratio.

For additional optimization tips, see [full protocol](#). Cell-type-specific recommendations available at [Reagent Agent: mirusbio.com/ra](http://mirusbio.com/ra)

# TransIT-X2<sup>®</sup> Dynamic Delivery System

## Quick Reference Protocol

*siRNA Delivery Instructions for MIR 6000, 6003, 6004, 6005, 6006, 6010*

*Full protocol, SDS and Certificate of Analysis available at [mirusbio.com/6000](http://mirusbio.com/6000)*



## SPECIFICATIONS

Storage	Store TransIT-X2 <sup>®</sup> Dynamic Delivery System tightly capped at -20°C. <b>Before each use</b> , warm to room temperature and vortex gently.
Product Guarantee	1 year from the date of purchase, when properly stored and handled.

### ► TRANSFECTION PROTOCOL | siRNA



Full protocol and additional documentation available at [mirusbio.com/6000](http://mirusbio.com/6000)

### Fill in volumes below based on culture vessel used for transfection (Table 2).

#### A. Plate cells

1. Plate cells in \_\_\_ml complete growth medium (per well or flask) 18-24 hours prior to transfection to ensure cells are actively dividing at the time of transfection.

**For adherent cells:** Plate cells at a density of 0.8—3.0 x 10<sup>5</sup> cells/ml.

**For suspension cells:** Plate cells at a density of 2.5—5.0 x 10<sup>5</sup> cells/ml.

2. Culture overnight. Most cell types should be ≥80% confluent at time of transfection.

#### B. Prepare TransIT-X2<sup>®</sup> Reagent:siRNA complexes

1. Warm TransIT-X2<sup>®</sup> Reagent to room temperature and vortex gently.
2. Place \_\_\_μl of OptiMEM<sup>®</sup> I Reduced-Serum Medium in a sterile tube.
3. Add \_\_\_μl TransIT-X2<sup>®</sup>. Mix gently by pipetting.
4. Add \_\_\_μl of a 10 μM siRNA stock solution (25 nM final concentration in well). Mix gently by pipetting.
5. Incubate at room temperature for 15-30 minutes.

#### C. Distribute complexes to cells

1. Add TransIT-X2<sup>®</sup>:siRNA complex mixture drop-wise to different areas of the well.
2. Gently rock plate for even distribution of complexes.
3. Incubate 24-72 hours.
4. Harvest cells and assay for knockdown of gene expression.

**Table 2.** Recommended starting conditions

Culture vessel	24-well plate	12-well plate	6-well plate
Surface area	1.9 cm <sup>2</sup>	3.8 cm <sup>2</sup>	9.6 cm <sup>2</sup>
Complete growth medium	0.5 ml	1 ml	2.5 ml
Serum-free medium	50 μl	100 μl	250 μl
siRNA (10 μM stock, 25 nM final)	1.4 μl	2.8 μl	6.8 μl
TransIT-X2 <sup>®</sup> Reagent	1.5 μl	3 μl	7.5 μl

### ► Transfection Optimization

Determine the best TransIT-X2<sup>®</sup>:DNA ratio for each cell type. Start with 3 μl of TransIT-X2<sup>®</sup> per 1 μg of DNA. Vary the concentration of TransIT-X2<sup>®</sup> from 2–6 μl per 1 μg DNA to find the optimal ratio.

*For additional optimization tips, see [full protocol](#). Cell-type-specific recommendations available at **Reagent Agent:** [mirusbio.com/ra](http://mirusbio.com/ra)*

# TransIT-X2<sup>®</sup> Dynamic Delivery System



## Quick Reference Protocol

DNA and siRNA Delivery Instructions for MIR 6000, 6003, 6004, 6005, 6006, 6010

Full protocol, SDS and Certificate of Analysis available at [mirusbio.com/6000](http://mirusbio.com/6000)

## SPECIFICATIONS

Storage	Store TransIT-X2 <sup>®</sup> Dynamic Delivery System tightly capped at -20°C. <b>Before each use</b> , warm to room temperature and vortex gently.
Product Guarantee	1 year from the date of purchase, when properly stored and handled.

### ▶ TRANSFECTION PROTOCOL | DNA & siRNA



Full protocol and additional documentation available at [mirusbio.com/6000](http://mirusbio.com/6000)

### Fill in volumes below based on culture vessel used for transfection (Table 3).

#### A. Plate cells

1. Plate cells in \_\_\_ ml complete growth medium (per well or flask) 18-24 hours prior to transfection to ensure cells are actively dividing at the time of transfection.

**For adherent cells:** Plate cells at a density of 0.8—3.0 x 10<sup>5</sup> cells/ml.

**For suspension cells:** Plate cells at a density of 2.5—5.0 x 10<sup>5</sup> cells/ml.

2. Culture overnight. Most cell types should be ≥80% confluent at time of transfection.

#### B. Prepare TransIT-X2<sup>®</sup> Reagent:DNA:siRNA complexes

1. Warm TransIT-X2<sup>®</sup> Reagent to room temperature and vortex gently.
2. Place \_\_\_ μl of OptiMEM<sup>®</sup> I Reduced-Serum Medium in a sterile tube.
3. Add \_\_\_ μl plasmid DNA. Mix gently by pipetting.
4. Add \_\_\_ μl of a 10 μM siRNA stock solution (25 nM final concentration in well). Mix gently by pipetting.
5. Add \_\_\_ μl TransIT-X2<sup>®</sup> Reagent. Mix gently by pipetting.
6. Incubate at room temperature for 15-30 minutes.

#### C. Distribute complexes to cells

1. Add co-transfection complex mixture drop-wise to different areas of the well.
2. Gently rock plate for even distribution of complexes.
3. Incubate 24-72 hours.
4. Harvest cells and assay as required.

**Table 3.** Recommended starting conditions

Culture vessel	24-well plate	12-well plate	6-well plate
Surface area	1.9 cm <sup>2</sup>	3.8 cm <sup>2</sup>	9.6 cm <sup>2</sup>
Complete growth medium	0.5 ml	1 ml	2.5 ml
Serum-free medium	50 μl	100 μl	250 μl
DNA (1 μg/μl stock)	0.5 μl	1 μl	2.5 μl
siRNA (10 μM stock, 25 nM final)	1.4 μl	2.8 μl	6.8 μl
TransIT-X2 <sup>®</sup> Reagent	1.5 μl	3 μl	7.5 μl

### ▶ Transfection Optimization

The amount of TransIT-X2<sup>®</sup> required for co-transfection is dictated by the amount of DNA. Determine the best TransIT-X2<sup>®</sup> Reagent:DNA ratio for each cell type. Start with 3 μl of TransIT-X2<sup>®</sup> per 1 μg of DNA. Vary the concentration of TransIT-X2<sup>®</sup> from 2–6 μl per 1 μg DNA to find the optimal ratio.

For additional optimization tips, see [full protocol](#). Cell-type-specific recommendations available at **Reagent Agent:** [mirusbio.com/ra](http://mirusbio.com/ra)

Mirus Bio LLC

[www.mirusbio.com](http://www.mirusbio.com) | [techsupport@mirusbio.com](mailto:techsupport@mirusbio.com) | U.S. Toll Free: 844.MIRUSBIO (844.647.8724) | Direct: +1.608.441.2852



**Reagent Agent<sup>®</sup>**

Reagent Agent<sup>®</sup> is an online tool designed to help determine the best solution for nucleic acid delivery based on in-house data, customer feedback and citations.

Learn more at: [mirusbio.com/ra](https://www.mirusbio.com/ra)

©1996-2024 All rights reserved. Mirus Bio LLC. All trademarks are the property of their respective owners.  
For terms and conditions, visit [www.mirusbio.com](https://www.mirusbio.com)

Rev.E 041019