

VirusGEN® AAV Transfection Kit

Quick Reference Protocol

Instructions for MIR 6745, 6750, 6755

Full protocol, SDS and Certificate of Analysis available at mirusbio.com/virusgen



SPECIFICATIONS

Storage	Store <i>TransIT</i> -VirusGEN® Transfection Reagent at -10 to -30°C, tightly capped. <i>Before each use</i> , warm to room temperature and vortex gently. Store VirusGEN® AAV Complex Formation Solution and Enhancer at 2 to 10°C.
Product Guarantee	When properly stored and handled, <i>TransIT</i> -VirusGEN® Transfection Reagent is guaranteed for 1 year from the date of purchase, and VirusGEN® Complex Formation Solution and Enhancer is guaranteed for 6 months from the date of purchase.

► PROTOCOL FOR ADENO-ASSOCIATED VIRUS (AAV) GENERATION IN SUSPENSION HEK 293 CULTURES



Full protocol and additional documentation available at mirusbio.com/virusgen

VirusGEN® AAV Transfection Kit Workflow

A. Maintain Cells
Passage cells regularly and ensure they are >95% viable before transfection.

B. Day 0
Dilute cells to $2 - 3 \times 10^6$ cells/ml.
Form transfection complexes in either VirusGEN® AAV CFS&E or PBS; use a volume that is 5-10% of cell culture volume.
Per ml of culture, add
1. **DNA:** 1 - 2 µg
2. **TransIT-VirusGEN®:** 1.5 - 3 µl.
Incubate stationary for **15 - 30 min**.

C. Add transfection complexes to cells.

D. Day 2 - 3
Harvest AAV **48 - 72 hr** post-transfection.

Fill in volumes below based on total culture volume (Table 1).

- A. Maintain cells**
1. Passage suspension HEK 293 cells 18-24 hours prior to transfection to obtain a density of $2 - 4 \times 10^6$ cells/ml the next day. Do NOT proceed with transfection if cells are not doubling every 24 hours or are < 95% viable.
 2. Incubate cells overnight at appropriate temperature and CO₂ levels.
- B. Prepare *TransIT*-VirusGEN® Reagent:DNA complexes**
1. At time of transfection, seed cells to a density of $2 - 3 \times 10^6$ cells/ml.
 2. Warm *TransIT*-VirusGEN® Reagent to room temperature and vortex.
 3. Place ____ml of VirusGEN® AAV Complex Formation Solution and Enhancer (CFS&E) or PBS in a sterile tube.
 4. Add ____µl of the total plasmid DNA to the tube. Mix gently by pipetting.
 5. Add ____µl of *TransIT*-VirusGEN® Reagent. Vortex gently to mix.
 6. Incubate at room temperature for 15-30 minutes to allow transfection complexes to form. Do not vigorously agitate complexes again once formed.
- C. Distribute complexes to cells in complete growth medium**
1. Add *TransIT*-VirusGEN® Reagent:DNA complexes (from Step B) to cells.
 2. Incubate cultures in appropriate conditions (i.e. 37°C, 5% CO₂, shaking) for 48-72 hours prior to AAV harvest.
- D. Harvest virus**
1. Following the 48-72 hour incubation, prepare 10X Cell Lysis Buffer (500 mM Tris pH 8, 10% Tween® 20, 20 mM MgCl₂).
 2. Transfer the total volume of cell suspension (____ml) to a sterile conical tube or appropriate vessel.
 3. Add 0.1X volume (____ml) of 10X Cell Lysis Buffer and 100 U/ml (____µl) of Benzonase®. Mix completely and incubate at 37°C for 1.5 hours with shaking.
 4. Add 0.1X volume (____ml) of 5 M NaCl. Mix completely and incubate at 37°C for 30 minutes with shaking.
 5. Centrifuge the mixture at 4,100 × g for 10 minutes to remove cell debris.
 6. Transfer the AAV-containing supernatant to a new tube. Store at -80°C.

Table 1. Volume scaling worksheet for VirusGEN® AAV Transfection Kit

Starting conditions per milliliter of complete growth medium				
	Per 1 ml	Total culture volume	Reagent quantities	
VirusGEN® AAV CFS&E or PBS	0.1 ml	× _____ml	=	_____ml
Total Plasmid DNA (1 µg/µl)	2 µl	× _____ml	=	_____µl
<i>TransIT</i> -VirusGEN® Reagent	3 µl	× _____ml	=	_____µl

Total Plasmid DNA refers to the combined mass of packaging plasmids and the transfer plasmid containing the gene-of-interest. Premix the plasmids together prior to adding to the complex formation medium.



Reagent Agent[®] is an online tool designed to help determine the best solution for nucleic acid delivery based on in-house data, customer feedback and citations.

Learn more at: mirusbio.com/ra

©1996-2025 All rights reserved. Mirus Bio LLC. All trademarks are the property of their respective owners.
For terms and conditions, visit www.mirusbio.com

Rev1 06142023

Mirus Bio LLC

www.mirusbio.com | techsupport@mirusbio.com | U.S. Toll Free: 844.647.8724 | Direct: +1.608.441.2852