

VirusGEN® LV Transfection Kit

Quick Reference Protocol

Instructions for MIR 6760, 6765, 6792

Full protocol, SDS and Certificate of Analysis available at mirusbio.com/virusgen




SPECIFICATIONS

Storage	Store <i>TransIT</i> -VirusGEN® Transfection Reagent at -10 to -30°C, tightly capped. Before each use , warm to room temperature and vortex gently. Store VirusGEN® LV Complex Formation Solution and VirusGEN® LV Enhancer at 2 to 10°C.
Product Guarantee	When properly stored and handled, <i>TransIT</i> -VirusGEN® Transfection Reagent is guaranteed for 1 year from the date of purchase, and VirusGEN® LV Complex Formation Solution and VirusGEN® LV Enhancer are guaranteed for 6 months from the date of purchase.

► PROTOCOL FOR LENTIVIRUS GENERATION IN ADHERENT HEK 293 CELL CULTURES



Full protocol and additional documentation available at mirusbio.com/virusgen



Fill in volumes below based on culture vessel used for transfection (Table 1).

A. Maintain cells

1. Approximately 18-24 hours prior to transfection, plate cells at a starting density of $4 - 6 \times 10^5$ cells/ml. Cultures should be 80-95% confluent at the time of transfection.
2. Incubate cells overnight at appropriate temperature and CO₂ levels.

B. Prepare *TransIT*-VirusGEN® Reagent:DNA complexes

1. Warm *TransIT*-VirusGEN® Reagent to room temperature and vortex.
2. Place ___ ml of VirusGEN® LV Complex Formation Solution in a sterile tube.
3. Add ___ μ l of the total plasmid DNA to the tube. Mix gently by pipetting.
4. Add ___ μ l of *TransIT*-VirusGEN® Reagent. Vortex gently to mix.
5. Incubate at room temperature for 15-30 minutes to allow transfection complexes to form. Do not vigorously agitate complexes again once formed.

C. Distribute complexes to cells in complete growth medium

1. Add *TransIT*-VirusGEN® Reagent:DNA complexes (from Step B) to cells.
2. Gently rock culture vessel for even distribution of complexes.
3. Incubate transfected cultures until addition of the VirusGEN® LV Enhancer.

D. Add VirusGEN® LV Enhancer to transfected cells

1. Between 5-24 hours post-transfection, add ___ mL of VirusGEN® LV Enhancer to cultures.
NOTE: Mirus recommends comparing transfection +/- Enhancer because, with some cell lines and vector constructs, equivalently high titer can be obtained *without* the Enhancer component.
2. Incubate cultures for an additional 24-43 hours prior to virus harvest (i.e. harvest lentivirus a total of 48 hours post-transfection).

E. Harvest virus

1. Following the 48-hour incubation, harvest cell supernatant containing recombinant lentivirus particles.
2. Filter through a 0.45 μ m PVDF filter to remove any cell debris.
3. Flash-freeze aliquots of lentivirus in cryotubes and store at -80°C.

Table 1. Recommended starting conditions

Culture Vessel	6-well plate	10-cm dish	T75 flask
Surface area	9.6 cm ²	59 cm ²	75 cm ²
Complete growth medium	2.0 ml	10 ml	15 ml
LV Complex Formation Solution	200 μ l	1.0 ml	1.5 ml
Total Plasmid DNA (1 μ g/ μ l)	2 μ l	10 μ l	15 μ l
<i>TransIT</i> -VirusGEN® Reagent	6 μ l	30 μ l	45 μ l
NOTE: Add VirusGEN® LV Enhancer 5-24 hours post-transfection.			
VirusGEN® LV Enhancer	200 μ l	1.0 ml	1.5 ml

Total Plasmid DNA refers to the combined mass of packaging plasmids and transfer plasmid. Premix the plasmids together prior to adding to complex formation medium.

Mirus Bio LLC

www.mirusbio.com | techsupport@mirusbio.com | U.S. Toll Free: 844.647.8724 | Direct: +1.608.441.2852

VirusGEN® LV Transfection Kit

Quick Reference Protocol

Instructions for MIR 6760, 6765, 6792

Full protocol, SDS and Certificate of Analysis available at mirusbio.com/virusgen



SPECIFICATIONS

Storage	Store <i>TransIT</i> -VirusGEN® Transfection Reagent at -10 to -30°C, tightly capped. Before each use , warm to room temperature and vortex gently. Store VirusGEN® LV Complex Formation Solution and VirusGEN® LV Enhancer at 2 to 10°C.
Product Guarantee	When properly stored and handled, <i>TransIT</i> -VirusGEN® Transfection Reagent is guaranteed for 1 year from the date of purchase, and VirusGEN® LV Complex Formation Solution and VirusGEN® LV Enhancer are guaranteed for 6 months from the date of purchase.

► PROTOCOL FOR LENTIVIRUS GENERATION IN SUSPENSION HEK 293 CELL CULTURES



Full protocol and additional documentation available at mirusbio.com/virusgen

Fill in volumes below based on total culture volume (Table 2).

A. Maintain cells

- Passage suspension HEK 293 cells 18-24 hours prior to transfection to obtain a density of $4 - 6 \times 10^5$ cells/ml the next day. Do NOT proceed with transfection if cells are not doubling every 24 hours or are < 95% viable.
- Incubate cells overnight (e.g. 37°C, 5-8% CO₂, shaking).

B. Prepare *TransIT*-VirusGEN® Reagent:DNA complexes

- Immediately prior to transfection, seed cells at a density of $2 - 4 \times 10^6$ cells/ml.
- Warm *TransIT*-VirusGEN® Reagent to room temperature and vortex.
- Place ___ ml of VirusGEN® LV Complex Formation Solution in a sterile tube.
- Add ___ µl of the total plasmid DNA to the tube. Mix completely.
- Add ___ µl of *TransIT*-VirusGEN® Reagent. Vortex gently to mix.
- Incubate for 15-30 minutes to allow transfection complexes to form. Do not vigorously agitate again once formed.

C. Distribute complexes to cells in complete growth medium

- Add *TransIT*-VirusGEN® Reagent:DNA complexes (from Step B) to cells.
- Incubate transfected cultures in an orbital shaker until addition of the VirusGEN® LV Enhancer.

D. Add VirusGEN® LV Enhancer to transfected cells

- Between 5-24 hours post-transfection, add ___ mL of VirusGEN® LV Enhancer to cultures.
- Incubate cultures for an additional 24-43 hours prior to virus harvest (i.e. harvest lentivirus a total of 48 hours post-transfection).

E. Harvest virus

- Following the 48-hour incubation, centrifuge cells in a sterile tube at $300 \times g$ for 5 minutes. Do NOT dispose of the supernatant.
- Transfer the virus-containing supernatant into a new sterile tube.
- Filter through a 0.45 µm PVDF filter to remove any cell debris.
- Flash-freeze aliquots of lentivirus in cryotubes and store at -80°C.

Table 2. Volume scaling worksheet for VirusGEN® LV Transfection Kit

Starting conditions per milliliter of complete growth medium			
	Per 1 ml	Total culture volume	Reagent quantities
LV Complex Formation Solution	0.1 ml	× _____ ml	= _____ ml
Total Plasmid DNA (1 µg/µl)	1 µl	× _____ ml	= _____ µl
<i>TransIT</i> -VirusGEN® Reagent	3 µl	× _____ ml	= _____ µl
NOTE: Add VirusGEN® LV Enhancer 18-24 hours post-transfection.			
VirusGEN® LV Enhancer	0.1 ml	× _____ ml	= _____ ml

Total Plasmid DNA refers to the combined mass of packaging plasmids and the transfer plasmid containing the gene-of-interest. Premix the plasmids together prior to adding to the complex formation medium.

©1996-2024 All rights reserved. Mirus Bio LLC. All trademarks are the property of their respective owners.

For terms and conditions, visit www.mirusbio.com

Rev1 060823

VirusGEN® LV Transfection Kit Workflow

A. Maintain Cells

Passage cells regularly; ensure >95% viability before transfection.

B. 0 hr

Dilute cells to $2 - 4 \times 10^6$ cells/ml.

Form complexes in VirusGEN® LV Complex Formation Solution in a volume that is 5-10% of cell culture volume.

Per ml of culture, add
1. DNA: 0.5 - 1.5 µg
2. *TransIT*-VirusGEN®: 1.5 - 4.5 µl.

Incubate stationary for **15 - 30 min.**

C. Add transfection complexes to cells.

D. 5 - 24 hr

Add **VirusGEN® LV Enhancer** in a volume that is 5-10% of cell culture volume.

E. 48 hr

Harvest lentivirus from supernatant.

Mirus Bio LLC

www.mirusbio.com | techsupport@mirusbio.com | U.S. Toll Free: 844.647.8724 | Direct: +1.608.441.2852